

Getting STARTED

Five steps to ensure your satisfaction:

1 DOWNLOAD "W&H BeHear" app




2 CONNECT BeHear to your phone

a. Make sure the smartphone Bluetooth connectivity is ON (Settings>>Bluetooth).


b. Press the BeHear NOW power button for about 5 seconds. LED will alternate red/blue.



For **iOS (iPhone) only:**
Go to iPhone "Settings" >> Bluetooth and add BeHear NOW to "My Devices".




- c. Launch the BeHear app.
- d. Search for **BeHear NOW** and tap the device name.
- e. The **BeHear NOW** should appear in **Paired** devices section.



3 HOW TO WEAR



Put the BeHear on your neck.
Insert earphones as shown.




Make sure "**W&H**" logo is right-side up.

4 CONDUCT hearing assessment

On app, tap **Personal** at the bottom of the screen.

Choose **Assessment**.

Follow the prompts, then press "**Apply**" and turn on the **BeHear** headset again.



5 CONTROL boxes

PERSONAL HEARING AMPLIFICATION
Right-hand buttons

- Increase ambient sounds
- Turn hearing mode ON/OFF (**long press**)
Switch hearing profiles (**short press**)
- Decrease ambient sounds



MUSIC PLAY/PHONE CALLS
Left-hand buttons

- Increase call/music volume
- Play/Pause audio (**short press**)
Accept calls/End calls (**short press**)
Reject calls (**long press**)
- Decrease call/music volume



NEED MORE ASSISTANCE? Let us help you!

Send us a message on our [Facebook Page](#)

E-mail us at support@wearandhear.com

Check out the FAQs on www.wearandhear.com

Need Assistance? E-mail us at support@wearandhear.com